

Anti-Human/Mouse/Rat ACBD6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02975

3th Edition

Description:ACBD6 (Acyl-CoA Binding Domain Containing 6) is a Protein Coding gene. Among its related pathways are Fatty Acyl-CoA Biosynthesis and Metabolism. GO annotations related to this gene include lipid binding and fatty-acyl-CoA binding. An important paralog of this gene is OSTF1.

Antigen:Synthesized peptide derived from the Internal region of human ACBD6

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Detected in placenta and spleen (at protein level). Detected in placenta, umbilical cord blood, CD34-positive hematopoietic progenitor cells and bone marrow.