

Anti-Human /Mouse RPA1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03008

3th Edition

Description:RPA1 (replication protein A1, 70kDa) is a protein-coding gene. Diseases associated with RPA1 include xeroderma pigmentosum, group a, and werner syndrome. GO annotations related to this gene include chromatin binding and single-stranded DNA binding.

Antigen:Recombinant protein of human RPA1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: