

**Anti-Human/Mouse/Rat Actin Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03017**

---

3th Edition

**Description:** Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

**Antigen:** Synthesized peptide derived from the C-terminal region of human Actin

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in some hepatocellular carcinomas