

Anti-Human/Mouse/Rat CORO1A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03024

3th Edition

Description: May be a crucial component of the cytoskeleton of highly motile cells, functioning both in the invagination of large pieces of plasma membrane, as well as in forming protrusions of the plasma membrane involved in cell locomotion. In mycobacteria-infected

Antigen: Synthesized peptide derived from the Internal region of human Coronin 1A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in brain, thymus, spleen, bone marrow and lymph node. Low in lung and gut.