

**Anti-Human /Mouse RNASE13 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03029**

---

3th Edition

**Description:** RNASE13 (Ribonuclease A Family Member 13 (Inactive)) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and ribonuclease activity.

**Antigen:** Recombinant protein of human RNASE13

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**