

Anti-Human/Mouse ADCK1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03039

3th Edition

Description:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).

Antigen:Synthesized peptide derived from the Internal region of human ADCK1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: