

Anti-Human COX2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03043

3th Edition

Description: Mediates the formation of prostaglandins from arachidonate. May have a role as a major mediator of inflammation and/or a role for prostanoid signaling in activity-dependent plasticity.

Antigen: Synthesized peptide derived from the C-terminal region of human Cox-2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: