

Anti-Human /Mouse /Rat EMD Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03058

3th Edition

Description:Emerin is a serine-rich nuclear membrane protein and a member of the nuclear lamina-associated protein family. It mediates membrane anchorage to the cytoskeleton. Dreifuss-Emery muscular dystrophy is an X-linked inherited degenerative myopathy resulting from mutation in the emerin gene.

Antigen:Recombinant protein of human EMD

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Skeletal muscle, heart, colon, testis, ovary and pancreas.