

Anti-Human AGBL4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03064

3th Edition

Description:AGBL4 (ATP/GTP Binding Protein Like 4) is a Protein Coding gene. Diseases associated with AGBL4 include Macular Degeneration, Age-Related, 1. GO annotations related to this gene include tubulin binding and metalloproteinase activity. An important paralog of this gene is AGBL5.

Antigen:Synthesized peptide derived from the C-terminal region of human AGBL4

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: