

Anti-Human/Mouse/Rat AK6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03084

3th Edition

Description: This gene encodes a protein that belongs to the adenylate kinase family of enzymes. The protein has a nuclear localization and contains Walker A (P-loop) and Walker B motifs and a metal-coordinating residue. The protein may be involved in regulation of Cajal body formation. In human, AK6 and TAF9 (GeneID: 6880) are two distinct genes that share 5' exons. Alternative splicing results in multiple transcript variants. AK6 (Adenylate Kinase 6) is a Protein Coding gene. Among its related pathways are Ribosome biogenesis in eukaryotes and Purine metabolism (REACTOME). GO annotations related to this gene include ATPase activity and adenylate kinase activity.

Antigen: Synthesized peptide derived from the N-terminal region of human AK6

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, chorionic villi and the central nervous system.