

**Anti-Human /Mouse /Rat HNRNPA2B1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03091**

---

3th Edition

**Description:**Involved with pre-mRNA processing. Forms complexes (ribonucleosomes) with at least 20 other different hnRNP and heterogeneous nuclear RNA in the nucleus.

**Antigen:**Recombinant protein of human HNRNPA2B1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**