

**Anti-Human /Mouse /Rat CD64 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03117**

---

3th Edition

**Description:** This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1.

**Antigen:** Recombinant protein of human FCGR1A

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Monocyte/macrophage specific.