

Anti-Human/Mouse/Rat CRMP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03121

3th Edition

Description:Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration.

Antigen:Synthesized peptide derived from human CRMP-2 around the non-phosphorylation site of T514.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous.