

Anti-Human/Mouse/Rat CstF-64T Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03174

3th Edition

Description: CSTF2T (Cleavage Stimulation Factor Subunit 2 Tau Variant) is a Protein Coding gene. Among its related pathways are mRNA Splicing - Major Pathway and Gene Expression. GO annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is CSTF2.

Antigen: Synthesized peptide derived from the Internal region of human CstF-64T

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: