

**Anti-Human/Mouse/Rat C-TAK1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03175**

---

3th Edition

**Description:**Involved in the specific phosphorylation of microtubule-associated proteins for tau, MAP2 and MAP4. Phosphorylates CDC25C on 'Ser-216'.

**Antigen:**Synthesized peptide derived from the N-terminal region of human C-TAK1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitous.