

**Anti-Human/Mouse/Rat/Monkey Cyclin A Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03198**

---

3th Edition

**Description:**Essential for the control of the cell cycle at the G1/S (start) and the G2/M (mitosis) transitions.

**Antigen:**Synthesized peptide derived from the Internal region of human Cyclin A

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Very high levels in testis and very low levels in brain. Also found in myeloid leukemia cell lines.