

Anti-Human/Monkey ALDH1B1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03234

3th Edition

Description: This protein belongs to the aldehyde dehydrogenases family of proteins. Aldehyde dehydrogenase is the second enzyme of the major oxidative pathway of alcohol metabolism. This gene does not contain introns in the coding sequence. The variation of this locus may affect the development of alcohol-related problems.

Antigen: Synthesized peptide derived from the Internal region of human ALDH1B1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Liver, testis and to a lesser extent in brain.