

Anti-Human/Mouse ALS2CR13 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03253

3th Edition

Description:FAM117B (Family With Sequence Similarity 117 Member B) is a Protein Coding gene. An important paralog of this gene is GLCC11.

Antigen:Synthesized peptide derived from the C-terminal region of human ALS2CR13

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: