

Anti-Human/Mouse/Rat ALX3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03258

3th Edition

Description: This gene encodes a nuclear protein with a homeobox DNA-binding domain that functions as a transcriptional regulator involved in cell-type differentiation and development. Preferential methylation of this gene's promoter is associated with advanced-stage neuroblastoma tumors. ALX3 (ALX Homeobox 3) is a Protein Coding gene. Diseases associated with ALX3 include Frontonasal Dysplasia 1 and Basal Encephalocele. GO annotations related to this gene include sequence-specific DNA binding. An important paralog of this gene is ALX4.

Antigen: Synthesized peptide derived from the Internal region of human ALX3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: