

Anti-Human/Mouse/Rat Amphiphysin I Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03269

3th Edition

Description: This gene encodes a protein associated with the cytoplasmic surface of synaptic vesicles. A subset of patients with stiff-man syndrome who were also affected by breast cancer are positive for autoantibodies against this protein. Alternate splicing of this gene results in two transcript variants encoding different isoforms. Additional splice variants have been described, but their full length sequences have not been determined. A pseudogene of this gene is found on chromosome 11.

Antigen: Synthesized peptide derived from the Internal region of human Amphiphysin I

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Neurons, certain endocrine cell types and spermatocytes.