

Anti-Human /Mouse /Rat CD83 Polyclonal Antibody
Polyclonal Antibody

Cat.NO.: PA03272

3th Edition

Description:The protein encoded by this gene is a single-pass type I membrane protein and member of the immunoglobulin superfamily of receptors. The encoded protein may be involved in the regulation of antigen presentation. A soluble form of this protein can bind to dendritic cells and inhibit their maturation. Three transcript variants encoding different isoforms have been found for this gene. CD83 (CD83 Molecule) is a Protein Coding gene. Diseases associated with CD83 include Exanthema Subitum and Specific Bursitis Often Of Occupational Origin. Among its related pathways are T cell receptor signaling pathway and B Cell Development Pathways.

Antigen:Recombinant protein of human CD83

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed by activated lymphocytes, Langerhans cells and interdigitating reticulum cells.