

Anti-Human/Mouse ANGPTL7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03277

3th Edition

Description:Angiotensin-converting enzyme 2 is a protein that in humans is encoded by the ANGPTL7 gene, a member of the ANGPTL family. Highly and specifically expressed in the cornea where is confined to the stromal layer.

Antigen:Synthesized peptide derived from the C-terminal region of human Angptl7

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highly and specifically expressed in the cornea where is confined to the stromal layer.