



Anti-Human NUSAP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03280

3th Edition

Description:NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization. May associate with chromosomes and promote the organization of mitotic spindle microtubules around them.

Antigen:Synthesized peptide derived from the C-terminal region of human ANKT

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: