

Anti-Human /Mouse /Rat LETMD1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03306

3th Edition

Description: This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene.

Antigen: Recombinant protein of human LETMD1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Kidney, liver, skeletal muscle, heart and brain. Overexpressed in various tumors including leukemia, lymphoma, and carcinomas of the breast, kidney, ovary, stomach, colon and uterine cervix.1 Publication