

**Anti-Human /Mouse /Rat 14-3-3 zeta/delta Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03345**

---

3th Edition

**Description:** Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

**Antigen:** Synthetic peptide of human YWHAZ(Ab-58)

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**