

**Anti-Human /Mouse FAM129B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03396**

---

3th Edition

**Description:**FAM129B (Family With Sequence Similarity 129 Member B) is a Protein Coding gene. Among its related pathways are MAP Kinase Signaling. An important paralog of this gene is FAM129A.

**Antigen:**Synthetic peptide of human FAM129B (Ab-712)

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**