

Anti-Human SLK Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03401

3th Edition

Description:SLK (STE20 Like Kinase) is a Protein Coding gene. Diseases associated with SLK include Amyotrophic Lateral Sclerosis 9. Among its related pathways are Apoptosis and Autophagy and Adhesion. GO annotations related to this gene include protein homodimerization activity and transferase activity, transferring phosphorus-containing groups. An important paralog of this gene is STK10.

Antigen:Synthetic peptide of human SLK

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed. Highest expression is found in heart and in skeletal muscle.