

**Anti-Human /Mouse /Rat MAP2K4 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03416**

---

3th Edition

**Description:**Dual specificity kinase that activates the JUN kinases MAPK8 (JNK1) and MAPK9 (JNK2) as well as MAPK14 (p38) but not MAPK1 (ERK2) or MAPK3 (ERK1).

**Antigen:**Synthetic peptide of human MAP2K4

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Abundant expression is seen in the skeletal muscle. It is also widely expressed in other tissues.