

Anti-Human /Mouse /Rat CD56 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03422

3th Edition

Description: This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.

Antigen: Synthetic peptide of human NCAM1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: