

Anti-Human /Mouse /Rat CD34 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03424

3th Edition

Description:The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene.

Antigen:Synthetic peptide of human CD34

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Selectively expressed on hematopoietic progenitor cells and the small vessel endothelium of a variety of tissues.