

Anti-Human/Monkey CYSLTR1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03433

3th Edition

Description: This gene encodes a member of the G-protein coupled receptor 1 family. The encoded protein is a receptor for cysteinyl leukotrienes, and is involved in mediating bronchoconstriction via activation of a phosphatidylinositol-calcium second messenger system. Activation of the encoded receptor results in contraction and proliferation of bronchial smooth muscle cells, eosinophil migration, and damage to the mucus layer in the lung. Upregulation of this gene is associated with asthma and dysregulation may also be implicated in cancer. Alternative splicing results in multiple transcript variants. CYSLTR1 (Cysteinyl Leukotriene Receptor 1) is a Protein Coding gene. Diseases associated with CYSLTR1 include Asthma and Adenoid Hypertrophy. Among its related pathways are RET signaling and Peptide ligand-binding receptors. GO annotations related to this gene include G-protein coupled receptor activity and leukotriene receptor activity. An important paralog of this gene is CYSLTR2.

Antigen: Synthesized peptide derived from the Internal region of human CysLTR1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed, with highest levels in spleen and peripheral blood leukocytes. Lower expression in several tissues, such as lung (mostly in smooth muscle bundles and alveolar macrophages), placenta, small intestine, pancreas, colon and heart.