

**Anti-Human CST11 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03436**

---

3th Edition

**Description:**Has antibacterial activity against the Gram-negative bacteria E.coli. May play a role in sperm maturation and fertilization.

**Antigen:**Synthesized peptide derived from the C-terminal region of human Cystatin 11

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Detected in the epithelium and lumen of the epididymis, and in sperm (at protein level).