

**Anti-Human/Mouse/Rat alpha Tubulin Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03448**

---

3th Edition

**Description:** Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

**Antigen:** Synthetic peptide of human TUBA1B

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**