

Anti-Human/Mouse/Rat YTHDF1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03491

3th Edition

Description:YTHDF1 (YTH N6-Methyladenosine RNA Binding Protein 1) is a Protein Coding gene. GO annotations related to this gene include poly(A) RNA binding and N6-methyladenosine-containing RNA binding. An important paralog of this gene is YTHDF3.

Antigen:Synthesized peptide derived from the N-terminal region of human DACA-1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: