

**Anti-Human/Mouse/Rat Desmin Polyclonal Antibody**  
**Polyclonal Antibody**

**Cat.NO.: PA03518**

---

3th Edition

**Description:** This gene encodes a muscle-specific class III intermediate filament. Homopolymers of this protein form a stable intracytoplasmic filamentous network connecting myofibrils to each other and to the plasma membrane. Mutations in this gene are associated with desmin-related myopathy, a familial cardiac and skeletal myopathy (CSM), and with distal myopathies.

**Antigen:** Synthesized peptide derived from human Desmin around the non-phosphorylation site of Ser60.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**