

Anti-Human OR51E1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03520

3th Edition

Description:Odorant receptor.

Antigen:Synthesized peptide derived from the C-terminal region of human D-GPCR

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highly expressed in prostate. Very low levels may be detected in some other tissues, such as placenta, skeletal muscle, heart, ovary and testis. Up-regulated in prostate cancers.