

**Anti-Human ZDHHC15 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03521**

---

3th Edition

**Description:**The protein encoded by this gene belongs to the DHHC palmitoyltransferase family. Mutations in this gene are associated with mental retardatio X-linked type 91 (MRX91). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. ZDHHC15 (Zinc Finger DHHC-Type Containing 15) is a Protein Coding gene. GO annotations related to this gene include receptor activity and palmitoyltransferase activity. An important paralog of this gene is ZDHHC2.

**Antigen:**Synthesized peptide derived from the C-terminal region of human DHHC-15

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in placenta, liver, lung, kidney, heart and brain.