

Anti-Human/Mouse/Rat DNAJB4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03551

3th Edition

Description: DNAJB4 belongs to the evolutionarily conserved DNAJ/HSP40 protein family. For background information on the DNAJ family, see MIM 608375.

Antigen: Synthesized peptide derived from the Internal region of human DnaJB4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in heart, pancreas and skeletal muscle, and to a lesser extent in brain, placenta and liver.