

Anti-Human/Mouse/Rat DSCR1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03608

3th Edition

Description: Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.

Antigen: Synthesized peptide derived from the C-terminal region of human DSCR 1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed heart, brain and skeletal muscle. Also expressed in all other tissues.