

Anti-Human DTYMK Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03612

3th Edition

Description:DTYMK (Deoxythymidylate Kinase) is a Protein Coding gene. Among its related pathways are superpathway of pyrimidine deoxyribonucleotides de novo biosynthesis and Purine metabolism (KEGG). GO annotations related to this gene include kinase activity and thymidylate kinase activity.

Antigen:Synthesized peptide derived from the N-terminal region of human dTMP Kinase

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: