

**Anti-Human DTYMK Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03612**

---

3th Edition

**Description:**DTYMK (Deoxythymidylate Kinase) is a Protein Coding gene. Among its related pathways are superpathway of pyrimidine deoxyribonucleotides de novo biosynthesis and Purine metabolism (KEGG). GO annotations related to this gene include kinase activity and thymidylate kinase activity.

**Antigen:**Synthesized peptide derived from the N-terminal region of human dTMP Kinase

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**