

Anti-Human/Mouse/Rat DCTN2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03620

3th Edition

Description: Modulates cytoplasmic dynein binding to an organelle, and plays a role in prometaphase chromosome alignment and spindle organization during mitosis. Involved in anchoring microtubules to centrosomes. May play a role in synapse formation during brain development.

Antigen: Synthesized peptide derived from the C-terminal region of human Dynactin 2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: