

**Anti-Human/Mouse/Rat EDG3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03648**

---

3th Edition

**Description:** Receptor for the lysosphingolipid sphingosine 1-phosphate (S1P). S1P is a bioactive lysophospholipid that elicits diverse physiological effect on most types of cells and tissues. When expressed in rat HTC4 hepatoma cells, is capable of mediating S1P-induced cell proliferation and suppression of apoptosis.

**Antigen:** Synthesized peptide derived from the Internal region of human EDG-3

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in all tissues, but most abundantly in heart, placenta, kidney, and liver.