

**Anti-Human/Mouse/Rat EDG5 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03651**

---

3th Edition

**Description:** This gene encodes a member of the G protein-coupled receptors, as well as the EDG family of proteins. This protein participates in sphingosine 1-phosphate-induced cell proliferation, survival, and transcriptional activation

**Antigen:** Synthesized peptide derived from the C-terminal region of human EDG-5

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**