

Anti-Human/Mouse/Rat EIF4E3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03711

3th Edition

Description: Recognizes and binds the 7-methylguanosine-containing mRNA cap during an early step in the initiation of protein synthesis and facilitates ribosome binding by inducing the unwinding of the mRNAs secondary structures.

Antigen: Synthesized peptide derived from the Internal region of human eIF4E3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: