

**Anti-Mouse/Rat PAFAH1B2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03721**

---

3th Edition

**Description:**Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene.

**Antigen:**Recombinant human Platelet activating factor acetylhydrolase 1b catalytic subunit 2 protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitous.