

**Anti-Mouse/Rat SIX1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03730**

---

3th Edition

**Description:**The protein encoded by this gene is a homeobox protein that is similar to the Drosophila 'sine oculis' gene product. This gene is found in a cluster of related genes on chromosome 14 and is thought to be involved in limb development. Defects in this gene are a cause of autosomal dominant deafness type 23 (DFNA23) and branchiootic syndrome type 3 (BOS3).

**Antigen:**Recombinant human Sine oculis homeobox homolog 1 protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Specifically expressed in skeletal muscle.