

Anti-Human/Mouse/Rat ELOVL6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03736

3th Edition

Description:Fatty acid elongases (EC 6.2.1.3), such as ELOVL6, use malonyl-CoA as a 2-carbon donor in the first and rate-limiting step of fatty acid elongation

Antigen:Synthesized peptide derived from the Internal region of human ELOVL6

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous.