

Anti-Mouse THRAP3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03740

3th Edition

Description:THRAP3 (Thyroid Hormone Receptor Associated Protein 3) is a Protein Coding gene. Among its related pathways are PEDF Induced Signaling and Circadian rythm related genes. GO annotations related to this gene include poly(A) RNA binding and transcription coactivator activity. An important paralog of this gene is BCLAF1.

Antigen:Recombinant human Thyroid hormone receptor-associated protein 3 protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous.