

**Anti-Human/Mouse/Rat/Monkey SH3GLB2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03750**

---

3th Edition

**Description:**SH3GLB2 (SH3 Domain Containing GRB2 Like, Endophilin B2) is a Protein Coding gene. Among its related pathways are Endocytosis. An important paralog of this gene is SH3GLB1.

**Antigen:**Synthesized peptide derived from the Internal region of human Endophilin B2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Detected in skeletal muscle, adipocyte, brain, lung, colon and mammary gland.